

ABSTRACT

Microarray technology provides the opportunity to identify thousands of microbial genes or populations simultaneously. Recently, a comprehensive functional gene array, called GeoChip 2.0, has been developed, evaluated and applied for characterizing microbial communities in natural systems. GeoChip 2.0 contains 24,243 oligonucleotide (50mer) probes and covers > 10,000 genes in >150 functional groups involved in nitrogen, carbon, sulfur and phosphorus cycling, metal reduction and resistance, and organic contaminant degradation. It is a powerful generic tool, and can be used for: (i) profiling various environmental samples, such as soil, groundwater, sediments, oil fields, deep sea, animal guts, and etc; (ii) studying biogeochemical processes and functional activities of microbial communities important to human health, agriculture, energy, global climate change, ecosystem management, and environmental cleanup and restoration; (iii) exploring direct linkages of microbial genes/populations to ecosystem processes and functions; and (iv) detecting functional genes and/or organisms in a particular environment. Here, we present an application example on the dynamics and stability of microbial genes and associated communities during a bioremediation period at the Oak Ridge Field Research Center (FRC). Due to exponential increases in the number of genes and the number of sequences for each gene, a new generation of such an array (GeoChip 3.0) is in development. GeoChip 3.0 is expected to have more features: (i) It is more comprehensive, covering >46,000 gene sequences of 292 gene families; (ii) it includes the phylogenetic marker, *gyrB*; (iii) It is more automatic for sequence retrieval and selection, probe design and verification, array construction and data analysis, information storage, and automatic update, which greatly facilitate the management of such a complicated array, especially for future updates.

Current version: GeoChip 2.0

Table 1 List of major functional markers on the GeoChip 2.0

Gene category	Example of key enzyme (gene)	Total probes
Nitrogen cycling		5309
Nitrogen fixation	Nitrogenase (<i>nifH</i>)	1225
Nitrification	Nitrite oxidoreductase (<i>nxrA</i>), nitrite oxidoreductase (<i>nxrB</i>), nitrite oxidoreductase (<i>nxrC</i>)	2366
Nitrification	Acetate monooxygenase (<i>amoA</i>), hydroxylamine monooxygenase (<i>hmoA</i>)	347
Nitrogen mineralization	Urease (<i>ureC</i>), glutamate dehydrogenase (<i>gdh</i>)	1432
Carbon cycling		4399
Carbon fixation	Rubisco (<i>rbcL</i> , <i>rbcS</i> , <i>rbcX</i>), C4H4, CO2H, PTHS	2018
Cellulose degradation	Cellulase, endoglucanase	1287
Lignin degradation	Ligninase, laccase	917
Chitin degradation	Chitinase, chitinase (<i>chiA</i>), chitinase	584
Methane production	Methyl coenzyme M reductase (<i>mcrA</i>)	437
Methane oxidation	Methanotrophs (<i>pmoA</i>)	136
Others	Lipase, peroxidase, lipase, peroxidase, cellulase	586
Sulfate reduction	Sulfate reductase (<i>srfA</i>), APS reductase	1613
Phosphorus utilization	Phosphatase (<i>phoA</i>)	145
Metal reduction and resistance		4546
Acetate resistance	Acetate reductase (<i>acrA</i> , <i>acrB</i>)	472
Chromium resistance	Chromium reductase (<i>crrA</i> , <i>crrB</i> , <i>crrC</i>)	262
Chromium resistance	Chromium reductase transporter (<i>crrA</i>)	139
Mercury resistance/reduction	Mercuric ion reductase (<i>merA</i> , <i>merB</i> , <i>merC</i>)	548
Nickel resistance	Nickel transporter (<i>nirA</i> , <i>nirB</i> , <i>nirC</i>)	140
Zinc resistance	Zinc reductase (<i>znrA</i>)	126
Other metal resistance	Other metal reductase, protein, sulfonamide reductase, etc.	2052
Contaminant degradation		8028
Benzene, toluene, chlorobenzene, and xylene (BTEX) & related aromatics	Benzene 1,2-dioxygenase (<i>bdoA</i>), chlorobenzene 1,2-dioxygenase (<i>cdoA</i>), toluene 1,2-dioxygenase (<i>tdoA</i>), xylene monooxygenase (<i>xmoA</i>), biphenyl CoA reductase (<i>bcrA</i>), and related 1,2-dioxygenase (<i>bdoA</i>)	4176
Chlorinated aromatics	Chlorobenzene reductase (<i>cbrA</i>)	95
Nitroaromatics	Nitrobenzene reductase (<i>nirA</i>), 4-nitrobenzoate reductase (<i>nirA</i>)	157
Polycyclic aromatic hydrocarbons (PAHs)	Polycyclic aromatic hydrocarbon reductase (<i>phrA</i>)	384
Polychlorinated biphenyls (PCBs)	Polychlorinated biphenyl reductase (<i>pcbA</i>)	761
Chlorinated alkenes (e.g. PCBs)	Chlorinated alkene reductase (<i>calA</i>)	252
Other organic compounds	Alkane hydroxylase (<i>alkA</i>), homocysteine 1,2-dioxygenase (<i>hdoA</i>), vanillin O-demethylase (<i>vdoA</i>)	2249
Total		24,243

CONCLUSIONS

- GeoChip has been constructed with more than 24,000 oligos covering more than 10,000 gene sequences. To our knowledge, this is the most comprehensive functional gene array currently available for environmental studies.
- GeoChip has been evaluated, and demonstrates that it can be used as a powerful tool for a rapid, high-throughput and cost-effective analysis of microbial communities.
- Microbial activities and associated communities were successfully monitored for *in situ* bioremediation at the Oak Ridge FRC site.
- A new generation of GeoChip (version 3.0) with more features is in development, which is expected to provide a more comprehensive picture for a given microbial community.

ACKNOWLEDGEMENT

ESPP2 (MDCase) is part of the Virtual Institute for Microbial Stress and Survival (VIMSS) supported by the U. S. Department of Energy, Office of Science, Office of Biological and Environmental Research, Genomics Program:GTL through contract DE-AC02-05CH11231 between Berkeley National Laboratory and the U. S. Department of Energy.

Further Development of GeoChip 3.0

New features for GeoChip 3.0

- GeoChip 3.0 is more comprehensive, and it contains >24,500 probes and covers about 47,000 (~10,000 for GeoChip 2.0) gene sequences of 292 gene families (~150 gene family on GeoChip 2.0) (Table 2). Thus, GeoChip 3.0 will be more representative.
 - The homology of automatically retrieved sequences by key words is verified by HUMMER using seed sequences so that unrelated sequences can be removed.
 - A software package (including databases) has been developed for sequence retrieval, probe and array design, probe verification, array construction, array data analysis, information storage, and automatic update, which greatly facilitate the management of such a complicated array, especially for future updates (Fig. 2).
 - GeoChip has implemented a universal standard, which can compare different samples, and normalize data.
 - GeoChip 3.0 implements a genomic control/standard, which can quantitatively analyze functional gene data.
- Automatic update greatly facilitates the management of such a complicated functional gene array.

Fig. 2 Work flow for GeoChip 3.0 design, construction and data analysis




Table 2. The summary of GeoChip 3.0 probe and sequence information

Gene category	No. of gene categories	No. of sequences for probe design	Total no. of probes designed	Total no. of CDS covered
Carbon degradation	31	9839	2720	4737
Carbon fixation	5	3378	898	1806
Methane reduction and oxidation	3	4182	254	434
Metal resistance and reduction	41	16825	4917	10458
Nitrogen cycling	13	27162	3561	6892
Organic remediation	190	31236	8815	16948
Phosphorus utilization	3	1441	599	1212
Sulfur cycling	3	4296	1328	1773
Energy process	2	901	413	449
Others (e.g. <i>gyrB</i>)	1	7957	1164	2251
Total	292	107217	24669	46960

Monitoring microbial activities during *in situ* bioremediation of uranium at the FRC site in Oak Ridge using GeoChip 2.0

Fig. 3 The microbial community dynamics from the monitoring well (102-2) during the bioremediation period from day 166 (3/1/2004) to day 719 (8/31/2005).




Fig. 4 Statistical analysis of richness across operational time and correlations (A) between GeoChip 2.0 results, geochemistry data, time and microbial communities (B, C, D).

